



■ NIH Chest X-ray14 Dataset

Exploratory Data Analysis Report

■ Dataset Overview

Metric	Value
Total Images	112,120
Unique Patients	30,805
Train Set	86,524 (77.2%)
Test Set	25,596 (22.8%)
Image Dimensions	~2518 × 2544 px
Multi-label Images	20,796 (18.5%)

■ Class Distribution



Disease Counts:


Class	Count	%
No Finding	60,361	53.8%

Infiltration	19,894	17.7%
Effusion	13,317	11.9%
Atelectasis	11,559	10.3%
Nodule	6,331	5.7%
Mass	5,782	5.2%
Pneumothorax	5,302	4.7%
Consolidation	4,667	4.2%
Pleural_Thick.	3,385	3.0%
Cardiomegaly	2,776	2.5%
Emphysema	2,516	2.2%
Edema	2,303	2.1%
Fibrosis	1,686	1.5%
Pneumonia	1,431	1.3%
Hernia	227	0.2%

■ ■ Class Imbalance Analysis


Imbalance Ratio: 266:1 (No Finding vs Hernia)

Recommendation: Use weighted loss or focal loss




■ ■ Multi-Label Analysis

- **Mean labels per image:** 1.26
- **Max labels:** 9 diseases in single image
- **Single-label:** 91,324 (81.5%)
- **Multi-label:** 20,796 (18.5%)







Disease Co-occurrence

Top combinations: Infiltration+Effusion (4,000), Effusion+Atelectasis (3,275)




■ Patient Demographics


- **Age:** Mean 46.9 years (range 1-95)
- **Gender:** Male 56.5%, Female 43.5%
- **View:** PA 60%, AP 40%



Disease Prevalence by Age




Disease Prevalence by Gender




Bounding Box Analysis

Warning: Only 984 bounding boxes (<1% coverage) - NOT suitable for YOLO!



Train/Test Split



■ Key Insights & Recommendations

Finding	Recommendation
Task Type	Multi-label Classification
Loss Function	Weighted BCE or Focal Loss
Output	Sigmoid (NOT Softmax)
Metric	AUC-ROC
Model	DenseNet-121 or ConvNeXt
YOLO?	■ NO - Insufficient BBox data

Suggested Models:

- 1. DenseNet-121 (CheXNet architecture - state-of-the-art)
- 2. ConvNeXt (modern architecture)
- 3. EfficientNet (good accuracy/efficiency)
- 4. Vision Transformers (ViT, Swin)